


Lecture 14: Linear Dimensionality Reduction for Spectra and Images

Instructor: Sergei V. Kalinin

Geometric Idea of PCA and LDA



Spectroscopic Imaging

Scanning probe microscopy:

- Force-distance curve measurements
- Current-voltage measurements
- Piezoresponse force/electrochemical strain spectroscopy

Electron microscopy:

- Electron Energy Loss Spectroscopy

Optical microscopy:

- Hyperspectral imaging
- Time resolved measurements

Mass-spectrometry:

- Secondary ion MS imaging

In many cases, measured signal can be represented or approximated as a linear combination of signals. However, their functional forms are generally unknown

Very important: convolution with resolution function is also mixing






Figure by
S. Jesse

Multiple Linear Regression

Linear mixing $S(\mathbf{x}, \mathbf{R}) = \sum_i a_i(\mathbf{x}) w_i(\mathbf{R}) + N$ but $w_i(\mathbf{R})$ are **known**


STEM of STO/LSMO/BFO interface Low-loss EELS spectra of three components



A.Y. BORISEVICH ET AL,
Suppression of Octahedral Tilts and Associated Changes in Electronic Properties at Epitaxial Oxide Heterostructure Interfaces, Phys. Rev. Lett. **105**, 087204 (2010).

“Chemistry”:
35 to 125 eV

Fit coefficient map





residuals map



χ^2 map



“Plasmons”
5 to 35 eV



Eigenvectors and loadings




Figure by A. Belianinov

Scree plot and correlations

- Semi log plot indicating the “weight” of each component as a function of all components
- Only the first few components contain useful info, while others are dominated by noise





Figure by A. Belianinov

- We can also analyze correlations in images

For AFM data

PCA Eigenvectors




Bayesian Linear Unmixing

$$S(\mathbf{x}, \mathbf{R}) = \sum_{i=1}^K a_i(\mathbf{x})w_i(\mathbf{R}) + N$$

$$\sum_{i=1}^K a_i(\mathbf{x}) = 1$$

- The eigenvectors $w_i(\mathbf{R})$ are non-negative, $w_i(\mathbf{R}) \geq 0$
- The loading coefficients sum to 1
- The number of eigenvectors, K , is a priori unknown




BLU is ideally suited for certain classes of problems, e.g. conduction through parallel channels, optical or electronic spectra of mixtures, etc


Laser heating induced phase transitions

- Copper indium thiophosphate ($\text{Cu}_{0.77}\text{In}_{1.12}\text{P}_2\text{S}_6$) layered ferroelectric
 - Ferroelectric state at room temperature
 - Curie temperature $T_c = 320$ K
 - Non-polar $\text{In}_{4/3}\text{P}_2\text{S}_6$ inclusions
- Combined Atomic Force Microscopy (AFM) and confocal Raman spectroscopy investigative approach
 - AFM – topography measurements
 - Piezoresponse force microscopy (PFM) – static ferroelectric domain structure
 - Raman – crystallographic structure via Raman spectra


Experimental scheme



Ferroelectric domain structure




Single point Raman spectra




A. IEVLEV, ACS Nano
9, 12442 (2015).

Laser heating induced phase transition


Laser can be used for local heating to induce ferroelectric- paraelectric phase transition



Domain structure evolution




Laser power



- Measurements with variation of the laser power
 - PFM – *in-situ* change in the domain structure above T_c
 - Raman – evolution of the Raman spectra through the phase transition
- Comprehensive analysis of Raman spectra is complicated by inhomogeneous chemical composition and high noise level
- Bayesian Linear Unmixing can be used for automated identification of spectra evolution


Raman spectra evolution (averaged)





BLU separation of components

Spatial concentration of components

Results of BLU: components and loading maps



A. IEVLEV, ACS Nano **9**, 12442 (2015).





Unmixing showed presence of three independent components in Raman spectra:


1. Non-polar $\text{In}_{4/3}\text{P}_2\text{S}_6$ – weak changes in intensity with temperature
2. Paraelectric CuInP_2S_6 above T_c – appears at higher laser powers
3. Ferroelectric CuInP_2S_6 below T_c – disappears at higher temperatures

PCA vs. LDA for elements

Alkali
Alkali
earth
Halogens





Alkali
Alkali
earth
Halogens





- Our element space is described by 13 descriptors
- In PCA, we found 2 linear combinations of these descriptors that describe this data set best.
- Alkali, alkali-earth, and halogens are close to each other in PCA space
- LDA finds best representation to separate alkali and halogens from everything else
- In LDA representation, alkali earth are close to alkali

What else can LDA give us?



- Decision function allows us to quantify how likely is the feature to belong to certain class

Visualizing the decision surfaces



- Generate multiple points uniformly distributed in the original high-dimensional space
- Perform the LDA transform
- Calculate the KDE